

OSG

List of EEG supplies for Braintronics




Table of Contents

| | |
|-------------------------------------|----|
| Braintronics accessories | 3 |
| Technomed electrodes / cables | 9 |
| EEG cup electrodes | 9 |
| Cables | 10 |
| Electrode linkers | 11 |
| Disposable EEG cup electrodes | 12 |
| Disposable EMG electrodes..... | 13 |
| Electro Caps and accessories..... | 14 |





Braintronics accessories



| Product code | Description |
|--------------|---|
| 0667/1BA | BrainBus termination connector for <u>EEG-11xx amplifiers</u>  |
| 0667/18BA | BrainNet termination connector for <u>EEG-1042 amplifier</u>  |
| 0667/2BA | Short BrainBus cable for connection of EEG-11xx to ISO-101  |
| 0667/6BA | 68-pin flat cable assy <u>15 cm</u> for EEG-11xx to INBOX-11xx  |
| 0667/7BA | 68-pin flat cable assy <u>1 mtr</u> for EEG-11xx to INBOX-11xx |
| 0667/8BA | BrainNet extension cable M8 female/male <u>2 mtr</u> (direct connection between headbox & FLASH) |

| | | |
|-----------|---|---|
| |  | |
| 0667/9BA | BrainNet extension cable M8 female/male <u>4 mtr</u> (direct connection between headbox & FLASH) |  |
| 0667/14BA | Short BrainNet cable with 5-pin large plug for easy detachment, ETHERNET-102 or ISO-101 side |  |
| 0667/19BA | BrainNet cable <u>6 mtr</u> with 5-pin large plug for easy detachment. (between two 14BA's) |  |
| 0667/20BA | BrainNet cable <u>10 mtr</u> with 5-pin large plug for easy detachment. (between two 14BA's) | |

| | | |
|---------------|---|---|
| |  | |
| 0667/21BA | PWR-24/15, Medical Grade adapter 24V / 15 Watt |  |
| 0667/21BA-IEC | Medical Grade adapter PWR-24/15 with short IEC male mains cable (voor aansluiting op isolatietransfo) |  |
| 0667/26BA | Event switch cable for EEG-1166, EEG-1142 or EEG-1042 amplifier |  |
| 0667/26BA-KR | Event switch cable for EEG-1166, EEG-1142 or EEG-1042 amplifier, curled cable | |


| | |
|-------------------|--|
| |  |
| 0667/26BA-KR-1.5M | Event switch cable for EEG-1166, EEG-1142 or EEG-1042 amplifier, curled cable 1.5M  |
| 0667/27BA | BrainNet cable 6 mtr for connection EEG-1042 to ETHERNET-102  |
| 0667/28BA | Block for fixation of EEG-1042 on microphone stand  |
| 0667/31BA | Adapter ElectroCap for EEG-1042, use cable 30BA or 32BA |

| | |
|-----------|--|
| |  |
| 0667/32BA | 40-pin flat cable assembly, short (approx. <u>15 cm</u>)  |
| 0667/30BA | 40-pin flat cable assembly. <u>1mtr.</u> |
| 0667/42BA | BrainBus Cable assembly for 2 Amplifier  |
| 0667/43BA | BrainBus Cable assembly for 3 Amplifier  |
| 0667/44BA | BrainBus Cable assembly for 4 Amplifier |

| | |
|-----------|---|
| |  <p>44BA BRAINBUS Cable Assembly for 4 Devices</p> <p>34BA Cable</p> <p>ISO-101 (V2)</p> |
| 0667/GND | <p>Ground Cable for Ethernet 102 Interface</p>  |
| 0667/4BB | <p>BrainBus Isolator for mini EEG amplifiers ISO-101</p>  |
| 0667/12BB | <p>Nonin X-pod with plug for EEG-1142/1166/1042 amplifier</p>  |

Technomed electrodes / cables

EEG cup electrodes



EEG cup electrodes have been designed to withstand long-term use by medical technicians. The precision stamped electrodes are available in two diameters: the pediatric, 6mm (0.24") or the adult, 10mm (0.39") size. They are stamped in two materials, 99.9% pure silver or tin. The silver electrodes are then covered with a heavy gold-plate, silver chlorided or kept silver. The flexible silicon coated wires are standard in multi colors.

| <u>CODE</u> | <u>PACKING</u> | <u>DESCRIPTION</u> |
|--------------|----------------------|--|
| 0234/C12-434 | per set of 12 pieces | (10 mm cup); 1 m silicon cable multi color with DIN 42 802 connector |
| 0234/C12-634 | per set of 12 pieces | (10 mm cup); 1,5 m silicon cable multi color with DIN 42 802 connector |
| 0234/C12-834 | per set of 12 pieces | (10 mm cup); 2 m silicon cable multi color with DIN 42 802 connector |
| 0234/C14-434 | per set of 12 pieces | (6 mm cup); 1 m silicon cable multi color with DIN 42 802 connector |
| 0234/C14-634 | per set of 12 pieces | (6 mm cup); 1,5 m silicon cable multi color with DIN 42 802 connector |
| <u>CODE</u> | <u>PACKING</u> | <u>DESCRIPTION</u> |
| 0234/C22-434 | per set of 12 pieces | (10 mm cup); 1 m silicon cable multi color with DIN 42 802 connector |
| 0234/C22-634 | per set of 12 pieces | (10 mm cup); 1,5 m silicon cable multi color with DIN 42 802 connector |
| 0234/C22-834 | per set of 12 pieces | (10 mm cup); 2 m silicon cable multi color with DIN 42 802 connector |

OSG B.V.B.A.
Bussestraat 17
B-2840 RUMST

Phone: +32 (0)15 32 13 73
Fax: +32 (0)15 32 13 93
Email: info@osg.be

| | | |
|--------------|----------------------|---|
| 0234/C24-434 | per set of 12 pieces | (6 mm cup); 1 m silicon cable multi color with DIN 42 802 connector |
| 0234/C24-634 | per set of 12 pieces | (6 mm cup); 1,5 m silicon cable multi color with DIN 42 802 connector |
| 0234/C24-834 | per set of 12 pieces | (6 mm cup); 2 m silicon cable multi color with DIN 42 802 connector |

| <u>CODE</u> | <u>PACKING</u> | <u>DESCRIPTION</u> |
|--------------|----------------------|--|
| 0234/C32-434 | per set of 12 pieces | (10 mm cup); 1 m silicon cable multi color with DIN 42 802 connector |
| 0234/C32-634 | per set of 12 pieces | (10 mm cup); 1,5 m silicon cable multi color with DIN 42 802 connector |
| 0234/C32-834 | per set of 12 pieces | (10 mm cup); 2 m silicon cable multi color with DIN 42 802 connector |
| 0234/C34-434 | per set of 12 pieces | (6 mm cup); 1 m silicon cable multi color with DIN 42 802 connector |
| 0234/C34-634 | per set of 12 pieces | (6 mm cup); 1,5 m silicon cable multi color with DIN 42 802 connector |
| 0234/C34-834 | per set of 12 pieces | (6 mm cup); 2 m silicon cable multi color with DIN 42 802 connector |

Cables

| Product code | Description |
|--------------|---|
| 0539/1041200 | <p>ECG cables: Reusable snap leads 5 colors, DIN 42 802 connector, 76 cm, per set of 5 pieces</p>  |

OSG B.V.B.A.
Bussestraat 17
B-2840 RUMST

Phone: +32 (0)15 32 13 73
Fax: +32 (0)15 32 13 93
Email: info@osg.be



| | |
|--------------|--|
| 0539/1075233 | <p>ECG cables: Reusable snap leads 5 colors, DIN 42 802 connector, 183 cm, per set of 5 pieces</p>  |
| 0539/1075400 | <p>EMG cables: Reusable snap leads, DIN 42 802 connector, 183 cm, per set of 10 pieces</p>  |
| 0539/1075401 | <p>EMG cables: Reusable snap leads, DIN 42 802 connector, 243 cm, per set of 10 pieces</p>  |
| 0234/D93-832 | <p>Reusable Extension cable 2m multicolor male-female DIN42802 connector 6 pc / package</p> |

Electrode linkers

| Product code | Description |
|--------------|-------------|
|--------------|-------------|

OSG B.V.B.A.
Bussestraat 17
B-2840 RUMST

Phone: +32 (0)15 32 13 73
Fax: +32 (0)15 32 13 93
Email: info@osg.be

| | |
|-----------------|--|
| 0234/ET30-101 | Electrode linker red |
| 0234/ET30-201 | Electrode linker black  |
| 0234/ET21-008RD | Double female adapter red, in: 2 mm pin connector, out: DIN 42 802 connector (per 25 pcs)  |
| 0234/ET22-008 | Female/male adapters, 2mm pin connector (per 25 pcs) |


Disposable EEG cup electrodes



| Product code | Packaging | Description |
|--------------|------------|--|
| 0234/C52-438 | 5x5 pieces | DISPOSABLE CUP AG/CL ELEKTRODES 10MM DIN42802 – 1meter |
| 0234/C52-838 | 5x5 pieces | DISPOSABLE CUP AG/CL ELEKTRODES 10MM DIN42802 – 2meter |





Disposable EMG electrodes





| Product code | Packaging | Description |
|--------------------|-----------|--|
| 0684/72010-K/10 | Per 420 | Neuroline 720 EMG elektrode, 22x30mm met kabel 100cm en K-connector (1,5mm) |
| 0684/72020-K/C/12 | Per 480 | Neuroline 720 EMG elektrode, 22x30mm met kabel 200cm en K-connector (1,5mm) |
| 0234/TE/K50430-002 | Per 40 | DISPOSABLE ADHESIVE SURFACE ELEKTRODES  |


Electro Caps and accessories

| Product code | Description |
|--------------|--|
| 0237/E1-L | Electro Cap Large (58-62 cm) blue  |
| 0237/E1-LM | Electro Cap L-M (56-60 cm) blue/red |
| 0237/E1-M | Electro Cap M (54-58 cm) red  |
| 0237/E1-MS | Electro Cap M-S (52-56 cm) red/yellow |
| 0237/E1-S | Electro Cap S (50-54 cm) yellow  |
| 0237/E1-SXS | Electro Cap S/XS (48-52 cm) yellow/green |
| 0237/E1-XS | Electro Cap XS (46-50 cm) green  |
| 0237/I 1-1 | Infra Cap I, 10/20 system, (42-46 cm) light blue |

| | |
|--------------|---|
| |  |
| 0237/I 1-2 | Infra Cap II, 10/10 system, (38-42 cm) pink  |
| 0237/I 1-2 V | Infra Cap II, 10/20 system, (38-42 cm) pink |
| 0237/I 1-3 | Infra Cap III, 10/10 system, (34-38 cm) brown  |
| 0237/I 1-3 V | Infra Cap III, 10/20 system, (34-38 cm) brown |
| 0237/I 1-4 | Infra Cap IV, 10/10 system, (30-34 cm) brown |
| 0237/I 1-5 | Infra Cap V, 10/10 system, (26-30 cm) brown |
| 0237/I 1-6 | Infra Cap VI, 10/10 system, (22-26 cm) brown |
| 0237/EC I-S | Electro Cap System I, with 1 Cap  <ul style="list-style-type: none">• Medium cap• Electrode board adaptor• Body harness medium• Quick insert electrode• Ear electrodes (pair)• Disposable sponge disks (pack)• Syringe / Needle kit• Electrode gel 16 Oz.• Special head measuring• Ivory liquid 12 Oz. |

| | |
|---------------|--|
| | <ul style="list-style-type: none"> • Instruction manual |
| 0237/EC II-S | <p>Electro Cap System II, with 2 Caps</p>  <ul style="list-style-type: none"> • Large cap • Medium cap • Electrode board adaptor • Body harness medium • 2x Quick insert electrode • 2x Ear electrodes (pair) • Disposable sponge disks (pack) • 2x Syringe / Needle kit • 2x Electrode gel 16 Oz. • Special head measuring • Ivory liquid 12 Oz. • Instruction manual |
| 0237/EC III-S | <p>Electro Cap System III, with 4 Caps</p>  <ul style="list-style-type: none"> • Large cap • Medium cap • Small cap • Extra small cap • Electrode board adaptor • Body harness medium • Body harness small / medium • 4x Quick insert electrode • 2x Ear electrodes (pair) • 2x Disposable sponge disks (pack) • 4x Syringe / Needle kit • 3x Electrode gel 16 Oz. • Special head measuring • Ivory liquid 12 Oz. |

| | |
|---------------|--|
| | <ul style="list-style-type: none"> • Instruction manual |
| 0237/INFA-SYS | Infra-Cap system with 3 Caps |
| 0237/E2S | Electrode board adapter, TP, 20 wires  |
| 0237/E4S | Quick insert electrode, TP  |
| 0237/E5-9S | Ear electrodes 9mm, pair, TP  |
| 0237/E5-6S | Ear electrodes 6mm, pair, TP  |
| 0237/E9 | Electro-gel 16 oz. (473ml)  |
| 0237/E13 | Cap straps (specify L, M, S or Extra Small)  |
| 0237/E23-5 | Chin straps 5" (12,5cm) |

| | |
|---------------|---|
| |  |
| 0237/E23-7 | Chin straps 7" (17,5cm) |
| 0237/10-30 | Nuprep 4 oz (114gr) pack of three |
| 0237/10-20-4 | Ten20 4 oz (114gr) jar, pack of three |
| 0237/10-20-4T | Ten20 4 oz (114gr) tube, pack of three |
| 0237/10-20-8 | Ten20 8 oz (228gr) pack of three |
| 0237/E7 | Needle(2)/Syringe(1) Kit,sterile |
| 0237/E6 | Disposable Sponge Disk package of 100 |