

OSG

---

# List of EEG supplies for Schwarzer AHNS



## Table of Contents




Schwarzer accessories.....	3
Sensors and cables.....	3
Connections on AHNS headbox .....	3
Technomed electrodes / cables .....	4
EEG cup electrodes .....	4
Cables .....	5
Electrode linkers .....	6
Disposable EEG cup electrodes .....	7
Disposable EMG electrodes.....	8
Electro Caps and accessories.....	9

# Schwarzer accessories

## Sensors and cables


### Connections on AHNS headbox



Product code	Description
0117/387018NE	Connection cable with Nonin Xpod Adapter 
0131/cap-AHNS	Cap connector for AHNS headbox 
0131/TEMP-409B	Temperature sensor: reusable attachable surface probe 

# Technomed electrodes / cables

## EEG cup electrodes



EEG cup electrodes have been designed to withstand long-term use by medical technicians. The precision stamped electrodes are available in two diameters: the pediatric, 6mm (0.24") or the adult, 10mm (0.39") size. They are stamped in two materials, 99.9% pure silver or tin. The silver electrodes are then covered with a heavy gold-plate, silver chlorided or kept silver. The flexible silicon coated wires are standard in multi colors.

<u>CODE</u>	<u>PACKING</u>	<u>DESCRIPTION</u>
0234/C12-434	per set of 12 pieces	(10 mm cup); 1 m silicon cable multi color with DIN 42 802 connector
0234/C12-634	per set of 12 pieces	(10 mm cup); 1,5 m silicon cable multi color with DIN 42 802 connector
0234/C12-834	per set of 12 pieces	(10 mm cup); 2 m silicon cable multi color with DIN 42 802 connector
0234/C14-434	per set of 12 pieces	(6 mm cup); 1 m silicon cable multi color with DIN 42 802 connector
0234/C14-634	per set of 12 pieces	(6 mm cup); 1,5 m silicon cable multi color with DIN 42 802 connector
<u>CODE</u>	<u>PACKING</u>	<u>DESCRIPTION</u>
0234/C22-434	per set of 12 pieces	(10 mm cup); 1 m silicon cable multi color with DIN 42 802 connector
0234/C22-634	per set of 12 pieces	(10 mm cup); 1,5 m silicon cable multi color with DIN 42 802 connector
0234/C22-834	per set of 12 pieces	(10 mm cup); 2 m silicon cable multi color with DIN 42 802 connector

OSG B.V.B.A.  
Bussestraat 17  
B-2840 RUMST

Phone: +32 (0)15 32 13 73  
Fax: +32 (0)15 32 13 93  
Email: [info@osg.be](mailto:info@osg.be)

0234/C24-434	per set of 12 pieces	(6 mm cup); 1 m silicon cable multi color with DIN 42 802 connector
0234/C24-634	per set of 12 pieces	(6 mm cup); 1,5 m silicon cable multi color with DIN 42 802 connector
0234/C24-834	per set of 12 pieces	(6 mm cup); 2 m silicon cable multi color with DIN 42 802 connector

<u>CODE</u>	<u>PACKING</u>	<u>DESCRIPTION</u>
0234/C32-434	per set of 12 pieces	(10 mm cup); 1 m silicon cable multi color with DIN 42 802 connector
0234/C32-634	per set of 12 pieces	(10 mm cup); 1,5 m silicon cable multi color with DIN 42 802 connector
0234/C32-834	per set of 12 pieces	(10 mm cup); 2 m silicon cable multi color with DIN 42 802 connector
0234/C34-434	per set of 12 pieces	(6 mm cup); 1 m silicon cable multi color with DIN 42 802 connector
0234/C34-634	per set of 12 pieces	(6 mm cup); 1,5 m silicon cable multi color with DIN 42 802 connector
0234/C34-834	per set of 12 pieces	(6 mm cup); 2 m silicon cable multi color with DIN 42 802 connector

## Cables

Product code	Description
0539/1041200	<p>ECG cables: Reusable snap leads 5 colors, DIN 42 802 connector, 76 cm, per set of 5 pieces</p> 

OSG B.V.B.A.  
Bussestraat 17  
B-2840 RUMST

Phone: +32 (0)15 32 13 73  
Fax: +32 (0)15 32 13 93  
Email: [info@osg.be](mailto:info@osg.be)



0539/1075233	<p>ECG cables: Reusable snap leads 5 colors, DIN 42 802 connector, 183 cm, per set of 5 pieces</p> 
0539/1075400	<p>EMG cables: Reusable snap leads, DIN 42 802 connector, 183 cm, per set of 10 pieces</p> 
0539/1075401	<p>EMG cables: Reusable snap leads, DIN 42 802 connector, 243 cm, per set of 10 pieces</p> 
0234/D93-832	<p>Reusable Extension cable 2m multicolor male-female DIN42802 connector 6 pc / package</p>

## Electrode linkers

Product code	Description
--------------	-------------

OSG B.V.B.A.  
Bussestraat 17  
B-2840 RUMST

Phone: +32 (0)15 32 13 73  
Fax: +32 (0)15 32 13 93  
Email: [info@osg.be](mailto:info@osg.be)

0234/ET30-101	Electrode linker red
0234/ET30-201	Electrode linker black 
0234/ET21-008RD	Double female adapter red, in: 2 mm pin connector, out: DIN 42 802 connector (per 25 pcs) 
0234/ET22-008	Female/male adapters, 2mm pin connector (per 25 pcs)


## Disposable EEG cup electrodes



Product code	Packaging	Description
0234/C52-438	5x5 pieces	DISPOSABLE CUP AG/CL ELEKTRODES 10MM DIN42802 – 1meter
0234/C52-838	5x5 pieces	DISPOSABLE CUP AG/CL ELEKTRODES 10MM DIN42802 – 2meter





## Disposable EMG electrodes





Product code	Packaging	Description
0684/72010-K/10	Per 420	Neuroline 720 EMG elektrode, 22x30mm met kabel 100cm en K-connector (1,5mm)
0684/72020-K/C/12	Per 480	Neuroline 720 EMG elektrode, 22x30mm met kabel 200cm en K-connector (1,5mm)
0234/TE/K50430-002	Per 40	DISPOSABLE ADHESIVE SURFACE ELEKTRODES 


## Electro Caps and accessories

Product code	Description
0237/E1-L	Electro Cap Large (58-62 cm) blue 
0237/E1-LM	Electro Cap L-M (56-60 cm) blue/red
0237/E1-M	Electro Cap M (54-58 cm) red 
0237/E1-MS	Electro Cap M-S (52-56 cm) red/yellow
0237/E1-S	Electro Cap S (50-54 cm) yellow 
0237/E1-SXS	Electro Cap S/XS (48-52 cm) yellow/green
0237/E1-XS	Electro Cap XS (46-50 cm) green 
0237/I 1-1	Infra Cap I, 10/20 system, (42-46 cm) light blue

	
0237/I 1-2	Infra Cap II, 10/10 system, (38-42 cm) pink 
0237/I 1-2 V	Infra Cap II, 10/20 system, (38-42 cm) pink
0237/I 1-3	Infra Cap III, 10/10 system, (34-38 cm) brown 
0237/I 1-3 V	Infra Cap III, 10/20 system, (34-38 cm) brown
0237/I 1-4	Infra Cap IV, 10/10 system, (30-34 cm) brown
0237/I 1-5	Infra Cap V, 10/10 system, (26-30 cm) brown
0237/I 1-6	Infra Cap VI, 10/10 system, (22-26 cm) brown
0237/EC I-S	Electro Cap System I, with 1 Cap  <ul style="list-style-type: none"><li>• Medium cap</li><li>• Electrode board adaptor</li><li>• Body harness medium</li><li>• Quick insert electrode</li><li>• Ear electrodes (pair)</li><li>• Disposable sponge disks (pack)</li><li>• Syringe / Needle kit</li><li>• Electrode gel 16 Oz.</li><li>• Special head measuring</li><li>• Ivory liquid 12 Oz.</li></ul>

	<ul style="list-style-type: none"> <li>• Instruction manual</li> </ul>
0237/EC II-S	<p>Electro Cap System II, with 2 Caps</p>  <ul style="list-style-type: none"> <li>• Large cap</li> <li>• Medium cap</li> <li>• Electrode board adaptor</li> <li>• Body harness medium</li> <li>• 2x Quick insert electrode</li> <li>• 2x Ear electrodes (pair)</li> <li>• Disposable sponge disks (pack)</li> <li>• 2x Syringe / Needle kit</li> <li>• 2x Electrode gel 16 Oz.</li> <li>• Special head measuring</li> <li>• Ivory liquid 12 Oz.</li> <li>• Instruction manual</li> </ul>
0237/EC III-S	<p>Electro Cap System III, with 4 Caps</p>  <ul style="list-style-type: none"> <li>• Large cap</li> <li>• Medium cap</li> <li>• Small cap</li> <li>• Extra small cap</li> <li>• Electrode board adaptor</li> <li>• Body harness medium</li> <li>• Body harness small / medium</li> <li>• 4x Quick insert electrode</li> <li>• 2x Ear electrodes (pair)</li> <li>• 2x Disposable sponge disks (pack)</li> <li>• 4x Syringe / Needle kit</li> <li>• 3x Electrode gel 16 Oz.</li> <li>• Special head measuring</li> <li>• Ivory liquid 12 Oz.</li> </ul>

	<ul style="list-style-type: none"> <li>• Instruction manual</li> </ul>
0237/INFA-SYS	Infra-Cap system with 3 Caps
0237/E2S	Electrode board adapter, TP, 20 wires 
0237/E4S	Quick insert electrode, TP 
0237/E5-9S	Ear electrodes 9mm, pair, TP 
0237/E5-6S	Ear electrodes 6mm, pair, TP 
0237/E9	Electro-gel 16 oz. (473ml) 
0237/E13	Cap straps (specify L, M, S or Extra Small) 
0237/E23-5	Chin straps 5" (12,5cm)

	
0237/E23-7	Chin straps 7" (17,5cm)
0237/10-30	Nuprep 4 oz (114gr) pack of three
0237/10-20-4	Ten20 4 oz (114gr) jar, pack of three
0237/10-20-4T	Ten20 4 oz (114gr) tube, pack of three
0237/10-20-8	Ten20 8 oz (228gr) pack of three
0237/E7	Needle(2)/Syringe(1) Kit,sterile
0237/E6	Disposable Sponge Disk package of 100