

OSG

---

# List of PSG supplies for Schwarzer AHNS



## Table of Contents

Schwarzer accessories.....	3
Sensors and cables.....	3
Use with PSG adapter (zwarte connector).....	3
Use with mini-headbox (safety connector).....	7
Connections on AHNS headbox .....	10
Velcro strap bands.....	11
Thermocan cannula/kabel.....	12
Nonin sensors.....	13
SleepSense inductive bands and sensors .....	15
Medcare locks .....	17
Technomed electrodes / cables .....	18
EEG cup electrodes .....	18
Cables.....	19
Electrode linkers .....	20
Disposable EEG cup electrodes .....	21
Disposable EMG electrodes.....	22
Pressure Transducer Airflow Sensor Pro-Tech.....	23
PTAF LITE.....	23
PTAF 2 .....	24
Cannulas (disposables) .....	24

## Schwarzer accessories





### Sensors and cables

Use with PSG adapter (zwarte connector)



Product code	Description
0117/593113	Snore sensor 
0117/593157 0539/P1758	Airflow sensor, 1 channel adult cannula style 
0117/593160	Airflow sensor for pediatrie 
0117/593153	Airflow sensor for babies

List of PSG supplies for Schwarzer AHNS

	
0117/503137	<p>Polygraphy adapter 20-channels</p> 
0117/635091	Cable PSG headbox - polygraphy adapter
0131/TCL-10	<p>Tibialis adapter lead with male DIN42802 connector, 10cm, 2 per package</p> <p>To be used with female DIN42802 – Ag/AgCl cup electrodes or female DIN42802 – Clip adapter lead</p> 
0117/555003	<p>2-poles EKG/ECG cable</p> 
0131/ECG-075	ECG connection lead, 75cm




		
0131/ECG-080	ECG connection lead, 80cm	
0131/ECG-100	ECG connection lead, 100cm	
0539/P1461	One Channel Kit Adult 1.5" Kit Contains: Adult Velcro® Strap Kit 1.5"; Piezo Sensor 1.5"  <b>(MET AANGEPASTE CONNECTOR)</b>	
		
0539/P1462	One Channel Kit Pediatric 1" Kit Contains: 1448: Pediatric VELCRO® Strap Kit 1"; 1453: Pediatric Piezo Sensor 1"  <b>(MET AANGEPASTE CONNECTOR)</b>	
		
0539/P1806	Pro-Tech zRIP Adult Effort Belt 	

List of PSG supplies for Schwarzer AHNS

0539/P1807	Pro-Tech zRIP Pediatric Effort Belt
0539/P1133	zRIP wireset <b>(MET AANGEPASTE CONNECTOR)</b>
0539/P1289	PTAF Lite Startup Kit <b>(met connection cable PTAF Lite - PSG adapter)</b> 
0131/PTAF-PSG	Connection cable PSG adapter - PTAF Lite – Flow sensor 

## Use with mini-headbox (safety connector)



Product code	Description
0539/P1716	Snore sensor 
0539/P1222	Airflow sensor 1 channel adult 
0539/P1223	Pro-Tech airflow sensor 1 channel pediatric 2 thermocouples nasal + 1 mouth 
0539/P1214	Pro-Tech airflow sensor, 1 channel pediatric small, nasal bridgemount

List of PSG supplies for Schwarzer AHNS


		
0539/P1273	Pro-Tech Airflow sensor 1 element	
0539/P1664 -WS2074	SPI DC body position sensor with safety connector	
0234/Z93-631	Reusable snap leads 5 colors, DIN 42 802 connector, 60" (150 cm), per set of 5 pieces	
0117/503140	Mini EB PSG/EEG 38 entrees, US Layout	



List of PSG supplies for Schwarzer AHNS



0117/503155	Mini EB EEG 38 entrees, US Layout
0117/503154	Mini EB EEG 38 entrees, Deutsche Layout
0117/635090	connection cable, mini-headbox – headbox (standaardlengte: 1,40m) (headbox met 3 rijen pinnen)
0117/635093	connection cable, mini-headbox – headbox (lengte: 3m) (headbox met 3 rijen pinnen)
0117/635092	adapter cable from 37pin headbox to 50pin mini-headbox (headbox met 2 rijen pinnen)
0539/P1461	One Channel Kit Adult 1.5" Kit Contains: Adult Velcro® Strap Kit 1.5"; Piezo Sensor 1.5" 
0539/P1462	One Channel Kit Pediatric 1" Kit Contains: 1448: Pediatric VELCRO® Strap Kit 1"; 1453: Pediatric Piezo Sensor 1" 
0539/P1806	Pro-Tech zRIP Adult Effort Belt 


## List of PSG supplies for Schwarzer AHNS

0539/P1807	Pro-Tech zRIP Pediatric Effort Belt
0539/P1133	zRIP wireset
0539/P1289	PTAF Lite Startup Kit 

## Connections on AHNS headbox





Product code	Description
0117/387018NE	Connection cable with Nonin Xpod Adapter 
0131/cap-AHNS	Cap connector for AHNS headbox 

0131/TEMP-409B	Temperature sensor: reusable attachable surface probe 
----------------	--

## Velcro strap bands

For use with:




- Old type of Schwarzer PSG adapter, ref. 210072
- New type of Schwarzer PSG adapter, ref. 503137
- Pro-Tech piezo sensor

Product code	Description
0117/860048 / 0539/P1420	Adult Velcro strap kit 1.5" (set = 3 belts) 
0539/P1421	Adult velcro strap, size 1
0539/P1422	Adult velcro strap, size 2
0539/P1423	Adult velcro strap, size 3
0117/860049 / 0539/P1448	Pediatric Velcro strap kit 1"(set = 4 belts) 


## Thermocan cannula/kabel

0454/15807	Thermocan 7-Ft Adult Thermal Cannula ( <u>per stuk</u> ) 
0454/1475	Adult Thermocan Interface Cable(Thermocouple) / Safety DIN 
0454/15813	Thermocan 7-Ft Pediatric Thermal Cannula ( <u>per stuk</u> ) 
0454/1475P	Pediatric Thermocan Interface Cable(Thermocouple) / Safety DIN 

## Nonin sensors

Product code	Description
0433/OX.8000J	<p>Nonin Adult Flex Sensor Includes 8000J sensor, attachment tape (25) and 1 meter cable</p> 
0433/OX.8000J-3	<p>Nonin Adult Flex Sensor Includes 8000J sensor, attachment tape (25) and <b>3 meter</b> cable</p>
0433/OX.8008J	<p>Nonin Infant Flex Sensor Includes 8008J sensor, Hydrogel tape strips (25) and 1 meter cable</p> 
0433/OX.8001J	<p>Nonin Neonatal Flex Sensor Includes 8001J sensor, Hydrogel tape strips (25) and 1 meter cable</p> 
0433/OX.8000JFW	Nonin flexi wrap adult attachment tape for Nonin 8000J sensor (25pcs)
0433/OX.8008JFW	Nonin flexi wrap pediatric attachment tape for Nonin 8008J sensor (25 pcs)
0433/OX.8000Q	Nonin ear clip sensor

# List of PSG supplies for Schwarzer AHNS

		
--	---	--

## SleepSense inductive bands and sensors



Product code	Description
0454/9002-L40	Multiple use inductive band - small (2/pack) 
0454/9002-L60	Multiple use inductive band - medium (2/pack) 
0454/9002-L90	Multiple use inductive band - large (2/pack) 
0454/9002-L120	Multiple use inductive band - x-large (2/pack) 
0454/9010	Reusable Inductive Effort Belt (Blue) ( <u>per piece</u> ) (1 year warranty) 
0454/9010-KIT	Reusable Inductive Effort Belt Kit (2 Belts + 9101 + 9102)
0454/9141	RIP Interface Cable 1.5 ft (45cm) - Abdomen / Key Connector ( <i>kabeltje voor aansluiting inductieve effort banden</i> ) 
0454/9142	RIP Interface Cable 1.5 ft (45cm) - Chest / Key Connector ( <i>kabeltje voor aansluiting inductieve effort banden</i> )

List of PSG supplies for Schwarzer AHNS

	
0454/9151	RIP Interface Cable <b>8 ft</b> (2,4m) - Abdomen / Key Connector
0454/9152	RIP Interface Cable <b>8 ft</b> (2.4m) - Chest / Key Connector
0454/9101	Inductive Interface Cable - Abdomen / Safety Connectors 
0454/9102	Inductive Interface Cable - Chest / Safety Connectors
0454/1250	Piezo Snore Sensor / Safety DIN connectors 
0454/14807	Luer to O2 adapter (pack of 5) 
0454/14808	Male to male Luer adapter(pack of 5) 




## Medcare locks

Product code	Description
0611/1421020_2	Medcare locks (XactTrace Lock – Abdomen + Thorax) 
0611/1421020	XactTrace Lock – Abdomen
0611/1421022	XactTrace Lock – Thorax
0611/1421000	XactTrace belt, 20m 

## Technomed electrodes / cables

### EEG cup electrodes



EEG cup electrodes have been designed to withstand long-term use by medical technicians. The precision stamped electrodes are available in two diameters: the pediatric, 6mm (0.24") or the adult, 10mm (0.39") size. They are stamped in two materials, 99.9% pure silver or tin. The silver electrodes are then covered with a heavy gold-plate, silver chlorided or kept silver. The flexible silicon coated wires are standard in multi colors.


<u>CODE</u>	<u>PACKING</u>	<u>DESCRIPTION</u>
0234/C12-434	per set of 12 pieces	(10 mm cup); 1 m silicon cable multi color with DIN 42 802 connector
0234/C12-634	per set of 12 pieces	(10 mm cup); 1,5 m silicon cable multi color with DIN 42 802 connector
0234/C12-834	per set of 12 pieces	(10 mm cup); 2 m silicon cable multi color with DIN 42 802 connector
0234/C14-434	per set of 12 pieces	(6 mm cup); 1 m silicon cable multi color with DIN 42 802 connector
0234/C14-634	per set of 12 pieces	(6 mm cup); 1,5 m silicon cable multi color with DIN 42 802 connector
<u>CODE</u>	<u>PACKING</u>	<u>DESCRIPTION</u>
0234/C22-434	per set of 12 pieces	(10 mm cup); 1 m silicon cable multi color with DIN 42 802 connector
0234/C22-634	per set of 12 pieces	(10 mm cup); 1,5 m silicon cable multi color with DIN 42 802 connector
0234/C22-834	per set of 12 pieces	(10 mm cup); 2 m silicon cable multi color with DIN 42 802 connector

## List of PSG supplies for Schwarzer AHNS

0234/C24-434	per set of 12 pieces	(6 mm cup); 1 m silicon cable multi color with DIN 42 802 connector
0234/C24-634	per set of 12 pieces	(6 mm cup); 1,5 m silicon cable multi color with DIN 42 802 connector
0234/C24-834	per set of 12 pieces	(6 mm cup); 2 m silicon cable multi color with DIN 42 802 connector

<u>CODE</u>	<u>PACKING</u>	<u>DESCRIPTION</u>
0234/C32-434	per set of 12 pieces	(10 mm cup); 1 m silicon cable multi color with DIN 42 802 connector
0234/C32-634	per set of 12 pieces	(10 mm cup); 1,5 m silicon cable multi color with DIN 42 802 connector
0234/C32-834	per set of 12 pieces	(10 mm cup); 2 m silicon cable multi color with DIN 42 802 connector
0234/C34-434	per set of 12 pieces	(6 mm cup); 1 m silicon cable multi color with DIN 42 802 connector
0234/C34-634	per set of 12 pieces	(6 mm cup); 1,5 m silicon cable multi color with DIN 42 802 connector
0234/C34-834	per set of 12 pieces	(6 mm cup); 2 m silicon cable multi color with DIN 42 802 connector

## Cables

<b>Product code</b>	<b>Description</b>
0539/1041200	ECG cables: Reusable snap leads 5 colors, DIN 42 802 connector, 76 cm, per set of 5 pieces 
0539/1075233	ECG cables: Reusable snap leads 5 colors, DIN 42 802 connector, 183 cm, per set of 5 pieces



## List of PSG supplies for Schwarzer AHNS

		
0539/1075400	EMG cables: Reusable snap leads, DIN 42 802 connector, 183 cm, per set of 10 pieces	
0539/1075401	EMG cables: Reusable snap leads, DIN 42 802 connector, 243 cm, per set of 10 pieces	
0234/D93-832	Reusable Extension cable 2m multicolor male-female DIN42802 connector 6 pc / package	

## Electrode linkers

Product code	Description
0234/ET30-101	Electrode linker red
0234/ET30-201	Electrode linker black

## List of PSG supplies for Schwarzer AHNS

	
0234/ET21-008RD	Double female adapter red, in: 2 mm pin connector, out: DIN 42 802 connector (per 25 pcs)
	
0234/ET22-008	Female/male adapters, 2mm pin connector (per 25 pcs)


## Disposable EEG cup electrodes



Product code	Packaging	Description
0234/C52-438	5x5 pieces	DISPOSABLE CUP AG/CL ELEKTRODES 10MM DIN42802 – 1meter
0234/C52-838	5x5 pieces	DISPOSABLE CUP AG/CL ELEKTRODES 10MM DIN42802 – 2meter

## Disposable EMG electrodes




Product code	Packaging	Description
0684/72010-K/10	Per 420	Neuroline 720 EMG elektrode, 22x30mm met kabel 100cm en K-connector (1,5mm)
0684/72020-K/C/12	Per 480	Neuroline 720 EMG elektrode, 22x30mm met kabel 200cm en K-connector (1,5mm)
0234/TE/K50430-002	Per 40	DISPOSABLE ADHESIVE SURFACE ELEKTRODES 

# Pressure Transducer Airflow Sensor Pro-Tech

## PTAF LITE

### SPECIFICATIONS:



- ✓ High and low level outputs
- ✓ Unfiltered airflow output
- ✓ Filtered airflow output
- ✓ Snoring output
- ✓ Self-test output signals
- ✓ Remote control options
- ✓ Auto and manual gain
- ✓ Auto Zero
- ✓ Differential pressure ports
- ✓ Filtered cannulas
- ✓ Battery life: 1000 hours

Product code	Description
0539/P1289	PTAF Lite Startup Kit

## PTAF 2

### SPECIFICATIONS:





- ✓ Unfiltered airflow output
- ✓ Filtered airflow output
- ✓ Snoring output
- ✓ Self-test output signals
- ✓ Remote control options
- ✓ Auto and manual gain
- ✓ Auto Zero
- ✓ Differential pressure ports
- ✓ Filtered cannulas
- ✓ Battery life: 1000 hours

Product code	Description
0539/P1291	PTAF2 Startup Kit (Requires a P1075 airflow cable and a P1075 snore cable) Kit Contains: -: 9V Alkaline Battery; P1258: Cannula Sample Pack; P1287: PTAF2 Sensor Module;
0539/P1075	Low Level Cable 1.5mm

## Cannulas (disposables)





Product code	Description
0539/P1259	ProTech adult sensor cannulas (60 per case)
0539/P1295	Pro-Tech adult nasal/oral cannulas (30 per case)
0539/P1263	Pro-Tech pediatric nasal sensor cannulas (60 per case)
0539/P1267	Pro-Tech pediatric small nasal sensor cannulas (60 per case)
0539/P1324	Cannula adapter kit (3 per set) <div data-bbox="459 734 743 904" data-label="Image">  </div>
0539/P1107	Cannula Titration Filter (50 pack) <div data-bbox="488 1021 654 1218" data-label="Image">  </div>