

OSG

List of PSG supplies






Table of Contents

Thermocan cannula/kabel.....	3
Nonin sensors.....	4
SleepSense inductive bands and sensors	6
Medcare locks	8
Technomed electrodes / cables	9
EEG cup electrodes	9
Cables	10
Electrode linkers	11
Disposable EEG cup electrodes	12
Disposable EMG electrodes.....	13
Pressure Transducer Airflow Sensor Pro-Tech	14
PTAF LITE.....	14
PTAF 2	15
Cannulas (disposables)	15
Schwarzer accessories.....	17
Sensors and cables.....	17
Use with PSG adapter (zwarte connector).....	17
Use with mini-headbox (safety connector).....	21
Connections on AHNS headbox	24
Velcro strap bands.....	25
Vitaport 4 PSG/CPS accessories	26
Morpheus accessories.....	29
Braintronics accessories	32


Thermocan cannula/kabel

0454/15807	Thermocan 7-Ft Adult Thermal Cannula (<u>per stuk</u>) 
0454/1475	Adult Thermocan Interface Cable(Thermocouple) / Safety DIN 
0454/15813	Thermocan 7-Ft Pediatric Thermal Cannula (<u>per stuk</u>) 
0454/1475P	Pediatric Thermocan Interface Cable(Thermocouple) / Safety DIN 

Nonin sensors

Product code	Description
0433/OX.8000J	<p>Nonin Adult Flex Sensor Includes 8000J sensor, attachment tape (25) and 1 meter cable</p> 
0433/OX.8000J-3	<p>Nonin Adult Flex Sensor Includes 8000J sensor, attachment tape (25) and 3 meter cable</p>
0433/OX.8008J	<p>Nonin Infant Flex Sensor Includes 8008J sensor, Hydrogel tape strips (25) and 1 meter cable</p> 
0433/OX.8001J	<p>Nonin Neonatal Flex Sensor Includes 8001J sensor, Hydrogel tape strips (25) and 1 meter cable</p> 
0433/OX.8000JFW	Nonin flexi wrap adult attachment tape for Nonin 8000J sensor (25pcs)
0433/OX.8008JFW	Nonin flexi wrap pediatric attachment tape for Nonin 8008J sensor (25 pcs)
0433/OX.8000Q	Nonin ear clip sensor





List of PSG supplies

		
--	---	--



SleepSense inductive bands and sensors

Product code	Description
0454/9002-L40	Multiple use inductive band - small (2/pack) 
0454/9002-L60	Multiple use inductive band - medium (2/pack) 
0454/9002-L90	Multiple use inductive band - large (2/pack) 
0454/9002-L120	Multiple use inductive band - x-large (2/pack) 
0454/9010	Reusable Inductive Effort Belt (Blue) (<u>per piece</u>) (1 year warranty) 
0454/9010-KIT	Reusable Inductive Effort Belt Kit (2 Belts + 9101 + 9102)
0454/9141	RIP Interface Cable 1.5 ft (45cm) - Abdomen / Key Connector (<i>kabeltje voor aansluiting inductieve effort banden</i>) 
0454/9142	RIP Interface Cable 1.5 ft (45cm) - Chest / Key Connector (<i>kabeltje voor aansluiting inductieve effort banden</i>)

List of PSG supplies


	
0454/9151	RIP Interface Cable 8 ft (2,4m) - Abdomen / Key Connector
0454/9152	RIP Interface Cable 8 ft (2.4m) - Chest / Key Connector
0454/9101	Inductive Interface Cable - Abdomen / Safety Connectors 
0454/9102	Inductive Interface Cable - Chest / Safety Connectors
0454/1250	Piezo Snore Sensor / Safety DIN connectors 
0454/14807	Luer to O2 adapter (pack of 5) 
0454/14808	Male to male Luer adapter(pack of 5) 

Medcare locks

Product code	Description
0611/1421020_2	Medcare locks (XactTrace Lock – Abdomen + Thorax) 
0611/1421020	XactTrace Lock – Abdomen
0611/1421022	XactTrace Lock – Thorax
0611/1421000	XactTrace belt, 20m 

Technomed electrodes / cables

EEG cup electrodes




EEG cup electrodes have been designed to withstand long-term use by medical technicians. The precision stamped electrodes are available in two diameters: the pediatric, 6mm (0.24") or the adult, 10mm (0.39") size. They are stamped in two materials, 99.9% pure silver or tin. The silver electrodes are then covered with a heavy gold-plate, silver chlorided or kept silver. The flexible silicon coated wires are standard in multi colors.

<u>CODE</u>	<u>PACKING</u>	<u>DESCRIPTION</u>
0234/C12-434	per set of 12 pieces	(10 mm cup); 1 m silicon cable multi color with DIN 42 802 connector
0234/C12-634	per set of 12 pieces	(10 mm cup); 1,5 m silicon cable multi color with DIN 42 802 connector
0234/C12-834	per set of 12 pieces	(10 mm cup); 2 m silicon cable multi color with DIN 42 802 connector
0234/C14-434	per set of 12 pieces	(6 mm cup); 1 m silicon cable multi color with DIN 42 802 connector
0234/C14-634	per set of 12 pieces	(6 mm cup); 1,5 m silicon cable multi color with DIN 42 802 connector
<u>CODE</u>	<u>PACKING</u>	<u>DESCRIPTION</u>
0234/C22-434	per set of 12 pieces	(10 mm cup); 1 m silicon cable multi color with DIN 42 802 connector
0234/C22-634	per set of 12 pieces	(10 mm cup); 1,5 m silicon cable multi color with DIN 42 802 connector
0234/C22-834	per set of 12 pieces	(10 mm cup); 2 m silicon cable multi color with DIN 42 802 connector

0234/C24-434	per set of 12 pieces	(6 mm cup); 1 m silicon cable multi color with DIN 42 802 connector
0234/C24-634	per set of 12 pieces	(6 mm cup); 1,5 m silicon cable multi color with DIN 42 802 connector
0234/C24-834	per set of 12 pieces	(6 mm cup); 2 m silicon cable multi color with DIN 42 802 connector

<u>CODE</u>	<u>PACKING</u>	<u>DESCRIPTION</u>
0234/C32-434	per set of 12 pieces	(10 mm cup); 1 m silicon cable multi color with DIN 42 802 connector
0234/C32-634	per set of 12 pieces	(10 mm cup); 1,5 m silicon cable multi color with DIN 42 802 connector
0234/C32-834	per set of 12 pieces	(10 mm cup); 2 m silicon cable multi color with DIN 42 802 connector
0234/C34-434	per set of 12 pieces	(6 mm cup); 1 m silicon cable multi color with DIN 42 802 connector
0234/C34-634	per set of 12 pieces	(6 mm cup); 1,5 m silicon cable multi color with DIN 42 802 connector
0234/C34-834	per set of 12 pieces	(6 mm cup); 2 m silicon cable multi color with DIN 42 802 connector



Cables

Product code	Description
0539/1041200	ECG cables: Reusable snap leads 5 colors, DIN 42 802 connector, 76 cm, per set of 5 pieces 
0539/1075233	ECG cables: Reusable snap leads 5 colors, DIN 42 802 connector, 183 cm, per set of 5 pieces

		
0539/1075400	EMG cables: Reusable snap leads, DIN 42 802 connector, 183 cm, per set of 10 pieces	
0539/1075401	EMG cables: Reusable snap leads, DIN 42 802 connector, 243 cm, per set of 10 pieces	
0234/D93-832	Reusable Extension cable 2m multicolor male-female DIN42802 connector 6 pc / package	

Electrode linkers

Product code	Description
0234/ET30-101	Electrode linker red
0234/ET30-201	Electrode linker black

	
0234/ET21-008RD	Double female adapter red, in: 2 mm pin connector, out: DIN 42 802 connector (per 25 pcs)
	
0234/ET22-008	Female/male adapters, 2mm pin connector (per 25 pcs)


Disposable EEG cup electrodes



Product code	Packaging	Description
0234/C52-438	5x5 pieces	DISPOSABLE CUP AG/CL ELEKTRODES 10MM DIN42802 – 1meter
0234/C52-838	5x5 pieces	DISPOSABLE CUP AG/CL ELEKTRODES 10MM DIN42802 – 2meter

Disposable EMG electrodes




Product code	Packaging	Description
0684/72010-K/10	Per 420	Neuroline 720 EMG elektrode, 22x30mm met kabel 100cm en K-connector (1,5mm)
0684/72020-K/C/12	Per 480	Neuroline 720 EMG elektrode, 22x30mm met kabel 200cm en K-connector (1,5mm)
0234/TE/K50430-002	Per 40	DISPOSABLE ADHESIVE SURFACE ELEKTRODES 

Pressure Transducer Airflow Sensor Pro-Tech

PTAF LITE

SPECIFICATIONS:



- ✓ High and low level outputs
- ✓ Unfiltered airflow output
- ✓ Filtered airflow output
- ✓ Snoring output
- ✓ Self-test output signals
- ✓ Remote control options
- ✓ Auto and manual gain
- ✓ Auto Zero
- ✓ Differential pressure ports
- ✓ Filtered cannulas
- ✓ Battery life: 1000 hours

Product code	Description
0539/P1289	PTAF Lite Startup Kit

PTAF 2

SPECIFICATIONS:





- ✓ Unfiltered airflow output
- ✓ Filtered airflow output
- ✓ Snoring output
- ✓ Self-test output signals
- ✓ Remote control options
- ✓ Auto and manual gain
- ✓ Auto Zero
- ✓ Differential pressure ports
- ✓ Filtered cannulas
- ✓ Battery life: 1000 hours

Product code	Description
0539/P1291	PTAF2 Startup Kit (Requires a P1075 airflow cable and a P1075 snore cable) Kit Contains: -: 9V Alkaline Battery; P1258: Cannula Sample Pack; P1287: PTAF2 Sensor Module;
0539/P1075	Low Level Cable 1.5mm

Cannulas (disposables)



Product code	Description
0539/P1259	ProTech adult sensor cannulas (60 per case)
0539/P1295	Pro-Tech adult nasal/oral cannulas (30 per case)
0539/P1263	Pro-Tech pediatric nasal sensor cannulas (60 per case)
0539/P1267	Pro-Tech pediatric small nasal sensor cannulas (60 per case)
0539/P1324	Cannula adapter kit (3 per set) 
0539/P1107	Cannula Titration Filter (50 pack) 





Schwarzer accessories




Sensors and cables

Use with PSG adapter (zwarte connector)



Product code	Description
0117/593113	Snore sensor 
0117/593157 0539/P1758	Airflow sensor, 1 channel adult cannula style 
0117/593160	Airflow sensor for pediatrie 
0117/593153	Airflow sensor for babies

	
0117/503137	<p>Polygraphy adapter 20-channels</p> 
0117/635091	Cable PSG headbox - polygraphy adapter
0131/TCL-10	<p>Tibialis adapter lead with male DIN42802 connector, 10cm, 2 per package</p> <p>To be used with female DIN42802 – Ag/AgCl cup electrodes or female DIN42802 – Clip adapter lead</p> 
0117/555003	<p>2-poles EKG/ECG cable</p> 
0131/ECG-075	ECG connection lead, 75cm




		
0131/ECG-080	ECG connection lead, 80cm	
0131/ECG-100	ECG connection lead, 100cm	
0539/P1461	One Channel Kit Adult 1.5" Kit Contains: Adult Velcro® Strap Kit 1.5"; Piezo Sensor 1.5" (MET AANGEPASTE CONNECTOR)	
		
0539/P1462	One Channel Kit Pediatric 1" Kit Contains: 1448: Pediatric VELCRO® Strap Kit 1"; 1453: Pediatric Piezo Sensor 1" (MET AANGEPASTE CONNECTOR)	
		
0539/P1806	Pro-Tech zRIP Adult Effort Belt 	

List of PSG supplies

0539/P1807	Pro-Tech zRIP Pediatric Effort Belt
0539/P1133	zRIP wireset (MET AANGEPASTE CONNECTOR)
0539/P1289	PTAF Lite Startup Kit (met connection cable PTAF Lite - PSG adapter) 
0131/PTAF-PSG	Connection cable PSG adapter - PTAF Lite – Flow sensor 

Use with mini-headbox (safety connector)




Product code	Description
0539/P1716	Snore sensor 
0539/P1222	Airflow sensor 1 channel adult 
0539/P1223	Pro-Tech airflow sensor 1 channel pediatric 2 thermocouples nasal + 1 mouth 
0539/P1214	Pro-Tech airflow sensor, 1 channel pediatric small, nasal bridgemount

		
0539/P1273	Pro-Tech Airflow sensor 1 element	
0539/P1664 -WS2074	SPI DC body position sensor with safety connector	
0234/Z93-631	Reusable snap leads 5 colors, DIN 42 802 connector, 60" (150 cm), per set of 5 pieces	
0117/503140	Mini EB PSG/EEG 38 entrees, US Layout	



0117/503155	Mini EB EEG 38 entrees, US Layout
0117/503154	Mini EB EEG 38 entrees, Deutsche Layout
0117/635090	connection cable, mini-headbox – headbox (standaardlengte: 1,40m) (headbox met 3 rijen pinnen)
0117/635093	connection cable, mini-headbox – headbox (lengte: 3m) (headbox met 3 rijen pinnen)
0117/635092	adapter cable from 37pin headbox to 50pin mini-headbox (headbox met 2 rijen pinnen)
0539/P1461	One Channel Kit Adult 1.5" Kit Contains: Adult Velcro® Strap Kit 1.5"; Piezo Sensor 1.5" 
0539/P1462	One Channel Kit Pediatric 1" Kit Contains: 1448: Pediatric VELCRO® Strap Kit 1"; 1453: Pediatric Piezo Sensor 1" 
0539/P1806	Pro-Tech zRIP Adult Effort Belt 


List of PSG supplies

0539/P1807	Pro-Tech zRIP Pediatric Effort Belt
0539/P1133	zRIP wireset
0539/P1289	PTAF Lite Startup Kit 

Connections on AHNS headbox





Product code	Description
0117/387018NE	Connection cable with Nonin Xpod Adapter 
0131/cap-AHNS	Cap connector for AHNS headbox 

0131/TEMP-409B	Temperature sensor: reusable attachable surface probe 
----------------	--


Velcro strap bands

For use with:

- Old type of Schwarzer PSG adapter, ref. 210072
- New type of Schwarzer PSG adapter, ref. 503137
- Pro-Tech piezo sensor







Product code	Description
0117/860048 / 0539/P1420	Adult Velcro strap kit 1.5" (set = 3 belts) 
0539/P1421	Adult velcro strap, size 1
0539/P1422	Adult velcro strap, size 2
0539/P1423	Adult velcro strap, size 3
0117/860049 / 0539/P1448	Pediatric Velcro strap kit 1"(set = 4 belts) 

Vitaport 4 PSG/CPS accessories





Product code	Description
0454/9010	Reusable Inductive Effort Belt (Blue) (<u>per piece</u>) (1 year warranty) 
0454/9002-L40	Multiple use inductive band - small (2/pack) 
0454/9002-L60	Multiple use inductive band - medium (2/pack) 
0454/9002-L90	Multiple use inductive band - large (2/pack) 
0454/9002-L120	Multiple use inductive band - x-large (2/pack) 
0131/SLP-patchpanel	Connection lead SLP – VP4 patchpanel (per set) (60 dagen garantie)

List of PSG supplies



		
0313/M92380-040	Airflow sensor, 1 channel neonatal nasal bridge mount, 40cm for VP4	
0313/M92385-040	Airflow sensor, 1 element	
0313/M92386-080	Airflow sensor, 1 channel adult cannula style (80cm cable)	
0313/M92387-040	Airflow sensor, 1 channel pediatric cannula style	
0313/M92506-040	Pro-Tech snore sensor for VP4 (40cm cable)	
0539/1041200-VP4	ECG snap lead set (3 leads x 76cm , red, yellow, green) for VP4 patchpanel	


List of PSG supplies

M98866.100	VP4 amplifier cable 1m 
M98866.020	VP4 amplifier cable 20cm
0313/M98830-G	Nonin Xpod VP4 ronde connector met single lock (1 gleufje), length 1m 
0313/M98830G.050	Nonin Xpod VP4 ronde connector met single lock (1 gleufje), length 50 cm






Morpheus accessories



Product code	Description
0602/ECG EL TPCLB	ECG Cable with connectors Touch Proof - Clip, length <u>150 cm</u> , Black 
0602/ECG EL TPCLR	ECG Cable with connectors Touch Proof - Clip, length <u>150 cm</u> , Red
0539/1041200-Morph	ECG cables: Reusable snap leads, DIN 42 802 connector, <u>76 cm</u> , per set of 2 pieces (black and white) 
0602/BSP OXY	Oxymeter to be used with MORPHEUS Recorder, with cable  
0602/BQ USB	USB Interface Card for MORPHEUS Serie Headbox 
0602/HBC INT LTM R	Cable to Connect Interfaces to MORPHEUS or to Quick Disconnect TERM LTM, New Model

	
0602/TERM LTM	Terminal to connect with MORPHEUS, with Quick Disconnect 
0602/BAG + BELT MORPH	Bag for Morpheus recorders + Belt for use with BAG MORPH  
0602/BQ STAND MORPH	MORPHEUS Headbox Stand, composed of BASE, STAND and Headbox Support SUPP MORPH 
0454/1250	Piezo Snore Sensor / Safety DIN connectors 
0539/P1222	Airflow sensor 1 channel adult

List of PSG supplies

		
0539/P1223	Pro-Tech airflow sensor 1 channel pediatric 2 thermocouples nasal + 1 mouth	
0539/P1214	Pro-Tech airflow sensor, 1 channel pediatric small, nasal bridgemount	
0539/P1273	Pro-Tech Airflow sensor 1 element	
0464/0514	Braebon cTherm Airflow Sensor (Cannula Compatible Airflow Sensor)	

Braintronics accessories

Product code	Description
0667/1BA	BrainBus termination connector for <u>EEG-11xx amplifiers</u> 
0667/18BA	BrainNet termination connector for <u>EEG-1042 amplifier</u> 
0667/2BA	Short BrainBus cable for connection of EEG-11xx to ISO-101 
0667/7BA	68-pin flat cable assy <u>1 mtr</u> for EEG-11xx to INBOX-11xx
0667/14BA	Short BrainNet cable with 5-pin large plug for easy detachment, ETHERNET-102 or ISO-101 side 
0667/19BA	BrainNet cable <u>6 mtr</u> with 5-pin large plug for easy detachment. (between two 14BA's)

		
0667/20BA	BrainNet cable <u>10 mtr</u> with 5-pin large plug for easy detachment. (between two 14BA's)	
0667/21BA	PWR-24/15, Medical Grade adapter 24V / 15 Watt	
0667/21BA-IEC	Medical Grade adapter PWR-24/15 with short IEC male mains cable (voor aansluiting op isolatietransfo)	
0667/26BA	Event switch cable for EEG-1166, EEG-1142 or EEG-1042 amplifier	
0667/27BA	BrainNet cable 6 mtr for connection EEG-1042 to ETHERNET-102	

List of PSG supplies

		
0667/28BA	Block for fixation of EEG-1042 on microphone stand 	
0667/GND	Ground Cable for Ethernet 102 Interface 	
0667/12BB	Nonin X-pod with plug for EEG-1142/1166/1042 amplifier 	